



CATHETER SELECTION & **CATHETERISATION**

Wilma Nicolson
BAUN Study Day

8th May 2015



PRINCIPLES OF CARE

- NHS Quality Improvement Scotland Best Practice Statement Urinary Catheterisation & Catheter Care Continence Adults with Urinary Dysfunction www.nhshealthquality.org
- RCN Catheter Care – www.rcn.org.uk
- Local NHS Urinary Catheterisation & Catheter Care Guidelines & Catheter Leaflets
- NES – Urinary Catheterisation, Management, Care & Prevention of Infection www.nes.scot.nhs.uk/hai/osc
- MHRA – Medicines and Healthcare products regulatory body – a medical device can be defined as any instrument, apparatus, appliance, material or healthcare product, excluding drugs, all catheter products require a CE mark.
- BS 1695 – Balloon size to fully inflate the balloon
- Report problems – Incident Forms Acute & Community
- Knowledge and Skills Framework
- Scottish Tariff www.isdscotland.org



PRINCIPLES OF CARE

- SIGN 88 – Management of Suspected Bacterial Urinary Tract Infections in Adults www.sign.ac.uk
- Health Protection Scotland Preventing Catheter Associated Urinary Tract Infection acute/community www.documents.hps.scot.nhs.uk
- European Association of Urology Nurses Catheterisation – Indwelling Catheters in Adults February 2012 – www.eaun.uroweb.org
- Local Acute & Hospital Catheter Formulary



NMC THE CODE 2015

- Maintain the knowledge and skills you need for safe and effective practice
- Keep to and promote recommended practice in relation to controlling infection
- Take all reasonable personal precautions necessary to avoid any potential health risks to colleagues, people receiving care and the public
- Record keeping – NMC



PRODUCT LIABILITY

- All catheters and equipment must be used within the manufacturers recommendations
- Any deviation from this will create liability on the part of the nurse





AIMS OF CATHETER CARE

- Promote patient comfort and dignity
- Recognise and minimise risk of CAUTI
- Assist patients achieve independence and self-care
- Able to provide cost effective service


TYPES OF URINARY CATHETERISATION

- Intermittent self catheterisation
- Indwelling catheterisation – catheter inserted into urethra
- Supra-pubic catheterisation – inserted into bladder surgically






CATHETERS

- A catheter is a thin hollow tube which can be inserted into the bladder through the urethra or supra-pubically
- Used by Romans, Chinese, Egyptians, material used reeds, gold, silver
- Allows for bacteria to enter the bladder
- Following insertion a biofilm forms on the surface of the catheter within hours and bacteria at the rate of 5% per day
- Within one month the bladder is fully colonised by bacteria
- Apart from trauma, infection accounts for virtually all complications of long term catheterisation



CATHETERS





CATHETER TYPES

- One way catheter – No balloon/intermittent – one channel for drainage, used for undertaking ISC or for treatment of urethral stricture
- Two way catheter – Balloon (Foley 1932) two internal channels, one for drainage and one to inflate the balloon (most commonly used)
- Three way balloon – Balloon (Foley) three channels – third to facilitate continuous bladder irrigation e.g. following urological surgery or if bleeding from bladder/prostate surgery



TYPES OF INTERMITTENT CATHETERISATION

- I.S.C. – Intermittent self catheterisation
- I.S.D. – Intermittent self dilation
- Via a cystostomy to drain urine



INTERMITTENT CATHETERISATION

- Insertion of a single use catheter into the bladder via the urethra to achieve regular and effective bladder emptying
- Preferred alternative to indwelling catheters when bladder emptying is incomplete – dependent on clinical assessment



INTERMITTENT CATHETERISATION

- Protects the kidneys, detrusor muscle & sphincter, prevents incontinence and decreases the number of urinary tract infections
- Improves quality of life and wellbeing – no indwelling catheter or legbag
- Allows patients to express their sexuality
- Should always be used in preference to indwelling catheter if clinically appropriate



REASON FOR INTERMITTENT SELF-CATHETERISATION

- Previous retention of urine
- MS
- CVA
- Parkinson's Disease
- Post Surgery



PATIENT CONSIDERATION

- Need to be prepared to do procedure
- Need compliance and ability
- Take equipment when you go out
- If not undertaking ISC can have increased incontinence/UTI – bladder overfills, back pressure, detrusor & sphincter damage
- Information available – DVD'S etc.
- Different types of catheters available



SELECTION OF INDWELLING CATHETER

- Reasons for catheterisation e.g. retention
- Patient consent, informed choice
- Individual needs of patient – who cares for catheter – patients/carer
- Sexually active/body image
- Allergy e.g. latex
- Length of time to remain in situ



SELECTION OF INDWELLING CATHETER

- Gender
- Size – smallest to allow drainage
- Comfort – easy to insert/remove – (silicone) prevent stricture
- Tissue compatibility
- **All catheters to be used within the manufacturers recommendations/information**

INDWELLING CATHETERS

– POINTS TO CONSIDER

- Material?
- Catheter size?
- Length of catheter?
- Balloon inflation volume?



TIP DESIGN



- Standard tip is round with two drainage eyes called a Nelaton catheter
- For routine catheterisation, a straight-tipped catheter should be used
- Tiemann has a covered tip, designed to negotiate the male prostatic curve, useful for difficult insertions
- Coude tip catheter has a curved tip like the Tiemann catheter, but has one, two or three drainage eyes situated in the curved tip



CATHETER MATERIAL

SHORT TERM P.T.F.E.



- Short/medium term – up to 28 days
- Refer to manufacturers instructions for urethral/supra-pubic use
 - Polytetrafluoroethylene Coated Latex (PTFE)
 - Do not use if any latex allergies
 - Smoother than plain latex which helps to prevent encrustation and irritation
 - Absorption of water is reduced due to the Teflon coating

CATHETER MATERIAL LONG TERM – HYDROGEL COATED

- Hydrogel coated latex – Up to 12 weeks
- Follow manufacturers guidelines for urethral or supra-pubic use
- Soft highly biocompatible
- Resistance to bacterial adherence
- Improved patient comfort
- Hydrophilic absorbs fluid to form a soft cushion around the catheter and reduces friction and urethral irritation, beneficial on insertion & removal



CATHETER MATERIAL LONG

TERM – SILICONE



- All silicone (100% silicone, hypoallergenic – uncoated, latex free) – up to 12 weeks
- Follow manufacturers guidelines for urethral or supra-pubic use
- Reduced tendency to encrustation
- While silicone causes less tissue damage, catheter has tendency to lose fluid, which increases risk of displacement in bladder
- Larger D shaped smooth wide lumen
- Greater risk of developing a cuff when deflated which can result in uncomfortable catheter removal or trauma


ANTIMICROBIAL & SILVER CATHETERS – COMMENTS

- Research has shown there is no clinical benefit in the use of these
- EUA state silver alloy and antibiotic impregnated catheters may decrease the risk of catheter associated bacteriurea in hospitalised patients
- No evidence that antibiotic-impregnated catheters decrease symptomatic infection and therefore they cannot be recommended routinely
- Cochrane review did not find sufficient evidence to determine the best type of catheter



CROSS SECTION OF CATHETER

- The inner lumen of the catheter varies quite a lot between catheter materials e.g. latex and a silicone catheter, so inserting a larger charriere catheter does not necessarily ensure a wider drainage channel





CATHETER SIZE

- Charriere – 19th century French instrument maker devised first sizing system
- One charriere unit is 1/3 of a mm
- A catheter marked 12ch has a 4.0mm external outer diameter
- Larger size more risk of urethral irritation and bypassing
- Catheters colour coded
- 6-10ch for paediatrics



CATHETER SIZE

- Routine drainage in an adult select smallest size which will effectively drain urine
- 12-14ch for normal urethral drainage – drain clear urine, no debris
- 14-18ch can be used for supra-pubic drainage – usually size 16fg or if debris/blood
- Size 20-24 – used for heavy haematuria, need for flushing (advice of urologist)
- Large sizes may be used for urine containing debris but will cause more urethral irritation and risk of bypassing



CATHETER LENGTH

- Paediatric length 30-31cm
- Consider sighting of legbag
- Female length 20-26cm (urethra 4-5cm) NOT TO BE USED FOR MALE CATHETERISATION
- Standard length – 40-45cm (male urethra 18-22cm) – also for supra-pubic, wheelchair/chair bound/overweight patient

INFLATION CHANNEL MARKINGS

- Product description
- Logo of company
- CE mark
- Balloon amount
- REF number
- Charriere (ch) size
- Outer diameter (mm)

CATHETER FUNNEL COLOR REFERENCE CHART

		5 fr		10 fr		16 fr		22 fr
		6 fr		12 fr		18 fr		24 fr
		8 fr		14 fr		20 fr		





CATHETER BALLOON SIZE

- **When catheter is in bladder, the balloon is inflated with sterile water/or glycerine to inflate balloon to keep catheter in bladder**
- **Balloon capacities are indicated by the inflation volume documented on the manufacturers instructions**
- **Some companies provide sterile pre-filled syringes with sterile water or glycerine 10%**
- **Fill with correct amount to stop distortion in bladder**
- **Under or over inflation can cause deflection of catheter tip and occlusion of eyelets of bladder mucosa, irritate the bladder wall and lead to bladder spasms**
- **Over inflation can cause stress on balloon**
- **Silicone catheters can lead to water loss from balloon over time**

CATHETER BALLOON SIZE

- For paediatric patients – 1.5ml, 3ml and 5ml
- 10ml – routine drainage (weighs 17gm)
- 30ml (weighs 48.2gm) – post prostatic surgery, can cause irritation of bladder lining, spasms, increased risk of infection due to drainage eyes sitting higher allowing for possible residual urine
- Use of larger balloon size will not prevent bypassing of urine
- Larger balloons tend to sit higher in the bladder with potential for increased residual urine volumes to collect below the catheter eyes



SELECTING A DRAINAGE SYSTEM/VALVE



■ Catheter Valves

- Catheter valves & drainage bags of all types should be considered to suit the patients needs
- Refer to local NHS guidelines
- Alternative to legbags
- Catheter valves are not always suitable for patients with detrusor instability, lack of bladder sensation or confused
- Small device which fits onto end of catheter
- Allows urine to be stored and drained from the bladder without the need for a permanently attached drainage bag
- Requires to be released at regular intervals to prevent over distension of bladder which may lead to UTIS or dilation of upper renal tract
- Can be used over 24 hours
- Can be attached to a night bag on free drainage



DRAINAGE BAGS

- **Assess patients dexterity, mobility, preference, cognitive function, lifestyle**
- **Ensure product concealed under clothing**
- **Who is going to empty/change leg/night bag**
- **Different materials/backings acute and community**
- **Must be emptied often enough to maintain urinary flow and prevent reflux**
- **A separate container must be used for each patient and contact between the tap and container avoided – follow local guidelines**



DAYTIME DRAINAGE BAGS

- **Anti reflux valve or anti-reflux chamber to prevent reflux of contaminated urine from the bag into the bladder**
- **Provide a closed drainage system, minimising risk of CAUTI**
- **Do not disconnect from catheter unless undertaking catheter maintenance solution or becomes contaminated**
- **Needle free sample port to obtain urine specimens**
- **Bag capacity can range from 150ml to 1000ml, the main capacity leg bags used are 350ml, 500ml, 750ml, empty according to patients daily routines**
- **Short/long tube this can vary from 5cm to 60cm, dependent on manufacturers**
- **Lever/t-tap/barrel**



DAYTIME DRAINAGE BAGS

- **Some companies pack contain pair of non latex gloves**
- **Straps come with leg bags**
- **Change every 7 days in line with manufacturers guidelines**
- **Never wash urine bags, do not add antiseptic or antimicrobial solutions to drainage bags**



NIGHT DRAINAGE BAGS

- **Single use 2 litres non-return twist off valve – connect to leg bag at night**
- **Drainable – if urine volumes require to be recorded or urinometer 2 litre or 3 litre**
- **Some companies e.g. LINC do a 2,600ml and Manfred Sauer 2,600mls**
- **Use of night stand – reduce risk of dislodging and improved drainage**



SPECIALIST DRAINAGE BAGS

- **Manfred Sauer**
 - Bendi bags, for wheelchair use
 - Discreet attached to thigh 450mls
 - Comfort range of bags, 600ml and 1000ml
 - 700-1330ml capacity
- **Rusch belly bags – body worn bag for urethral and supra-pubic use or nephrostomy, has an anti-reflux valve behind the catheter port which prevents reflux urine flow which allows positioning of this bag above level of bladder**





CATHETER SECUREMENT DEVICES

- **Prevents excess traction of catheter against bladder neck**
- **Securement of the catheter reduces catheter movement, infection risk and improves patient comfort**



NON-ADHESIVE DEVICES

- **Leg straps – follow manufacturers instructions**
- **Sleeves – measure thigh circumference**
- **G-straps/Comfrasure – follow manufacturers instructions**



CATHETER SECUREMENT – ADHESIVE DEVICES



- **Statlock – adhesive foam pad with swivel clip and locking device**
- **Clinifix – hydrocolloid strip – can be fixed or allows some movement**
- **Follow manufacturers guidelines**
- **On prescription or PECOS**






DOCUMENTATION

- **Peelable record labels that record catheter material, manufacturer, lot/batch number, size of catheter**
- **Amount of water in balloon**
- **Problem with insertion/removal**
- **Reason for catheterisation/re-catheterisation**
- **Instilligel used amount, expiry date, lot/batch number**
- **Catheter information leaflet give**
- **Signature of nurse**





Any Questions

The End